

BlueArctic

cryotherapy

continuous cold

Key Features and Benefits

- Two position locking top with carrying handle
- Water tight seal
- Hose stirrup
- Inlet temperature control & monitoring
- Semi-closed loop system
- Durable and dependable
- Compatibility with a variety of pads

Temperature Performance

The ideal temperature range for cold therapy is 40° - 50°F. The BlueArctic maintains a near-constant temperature in this range for up to 7 hours of comfort and therapy for your patients. With it's patented semi-closed loop system, the Blue Arctic achieves a more constant and accurate temperature than the competitor's systems on a single charge of ice.

Universal Quick Connectors

The Blue Arctic's easy pad connectors allow patients to be quickly mobile. They also allow for a wide variety of of cold therapy pads. OrthoRehab can order sterile pads to meet any size and application need you may have.



The Blue Arctic Cold Therapy unit helps reduce the pain, haemarthrosis, edema and inflammation which hinders rehabilitation and compliance.

passionate about motion.™

ORTHO REHAB™



BlueArctic

cryotherapy

continuous cold

45°

Why Use Continuous Cold Therapy?

Cooling not only decreases nerve-ending excitability, it also slows metabolism and reduces muscle spasms. The amount of time required to achieve significant muscle cooling can vary depending upon the thickness of the subcutaneous fat layer. The longer the exposure, the greater the depth of cooling effect. However, direct ice application for non-controlled cold therapy can cause adverse reactions such as frostbite and nerve palsy when the temperature is below the safe range and/or applied for too long.

With the safe, controlled temperature of the BlueArctic your patients have hours of low-maintenance therapy that can help to reduce bleeding, edema, inflammation and the need for pain medications.

Contraindications:

Cold therapy should not be used on persons with:

- Raynaud's or other vasospastic diseases
- Sickle-Cell Anemia
- Cryoglobulinemia
- Cold hypersensitivity
- Compromised local circulation.

Discontinue use and consult your physician immediately if numbness or skin reddening occurs.

Ordering Information:

PART NUMBER

BlueArc BlueArctic Continuous Cold Therapy Unit

0693-9 Battery Pack

0674-9 Dual Pad Connector

0612-9 Cold Wrap Shoulder, Sterile L/XL

0801-9 Cold Pad Universal, Sterile

Ideal Application

The BlueArctic Cold Therapy unit helps reduce pain and swelling while facilitating rehabilitation. It provides up to 7 hours of continuous cold therapy in the ideal temperature range for the maximum benefit to patients.

ORTHOREHAB™

In the USA call:

1.877.4.MOTION

(-877-466-8466) to contact
your local sales or service representative

www.orthorehab.com

International Customers contact:

OrthoMotion Inc.

901 Dillingham Road

Pickering, Ontario L1W 2Y5 Canada

1.905.420.3303;

E-mail: customerservice@orthomotion.ca

For more information contact your local representative